

Anti-CD210 Rat Monoclonal Antibody, FITC-Cy5 Conjugated



Catalog #: F0770-FITC

Aliases

IL10RA; Interleukin 10 Receptor Subunit Alpha; HIL-10R; CDW210A; CD210a; CD210; IL10R; Interleukin-10 Receptor Subunit Alpha; Interleukin-10 Receptor Subunit 1; Interleukin 10 Receptor, Alpha; IL-10 Receptor Subunit Alpha; IL-10R Subunit Alpha; IL-10R Subunit 1; IL-10R1; IL-10RA; Interleukin-10 Receptor Alpha Chain; CD210 Antigen; CDw210a

Background

Gene Name: IL10RA

NCBI Gene Entry: [3587](#)

UniProt Entry: [Q13651](#)

Application Information

Clonality: Rat monoclonal antibody

Clone ID: 24GB15685

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Rat IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD210 Rat Monoclonal Antibody, FITC- Cy5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM