

Anti-FOXP3 Mouse Monoclonal Antibody, FITC-Cy5 Conjugated



Catalog #: F0772-FITC

Aliases

FOXP3; Forkhead Box P3; DIETER; XPID; AIID; PIDX; IPEX; JM2; Immune Dysregulation, Polyendocrinopathy, Enteropathy, X-Linked; Forkhead Box Protein P3; SCURFIN; Immunodeficiency, Polyendocrinopathy, Enteropathy, X-Linked; FOXP3delta7; Scurfin

Background

Gene Name: FOXP3

NCBI Gene Entry: [50943](#)

UniProt Entry: [Q9BZS1](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15695

Species Reactivity: Mouse

Applications Tested: Intracellular Staining for Flow Cytometry ☐ ICFC ☐

Conjugation: FITC-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG1, κ

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-FOXP3 Mouse Monoclonal Antibody, FITC-Cy5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM