

Anti-H-2Kd/H-2Dd Mouse Monoclonal Antibody, FITC-Cy5 Conjugated



Catalog #: F0775-FITC

Aliases

H-2D; H2-D; K-f; H-2K; H2-K; H2-D1; H-2K(d)

Background

Gene Name: H2-D1/H2-K1

NCBI Gene Entry: [14964/14972](#)

Application Information

Clonality: Mouse monoclonal antibody

Clone ID: 24GB15700

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Mouse IgG2a

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-H-2Kd/H-2Dd Mouse Monoclonal Antibody, FITC-Cy5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM