

Anti-Ly-6G/Ly-6C Rat Monoclonal Antibody, FITC-Cy5 Conjugated



Catalog #: F0780-FITC

Aliases

Gr1; Gr-1; Ly-6G; Ly6c; Ly-6C; Ly-6C1

Background

Gene Name: Ly6g/Ly6c1

NCBI Gene Entry: [546644/17067](#)

Application Information

Clonality: Rat monoclonal antibody

Clone ID: 24GB15630

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Rat IgG2b

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Ly-6G/Ly-6C Rat Monoclonal Antibody, FITC-Cy5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM