

Anti-TCR β Armenian Hamster Monoclonal Antibody, FITC-Cy5 Conjugated



Catalog #: F0783-FITC

Aliases

TRB; T Cell Receptor Beta Locus; T-Cell Antigen Receptor, Beta Polypeptide, T-Cell Receptor, Beta Cluster; M1-Specific T Cell Receptor Beta Chain; TR Beta Chain
TRBV19*01J2S7*01C*02; TCRB, TRB@; TR Beta Chain TRBV25-1*01J2S3*01C2*01; T Cell Receptor Beta Chain MC.7.G5; T-Cell Receptor, Beta Cluster; MC.7.G5 TRB; TCRB; TRB@

Background

Gene Name: TRB

NCBI Gene Entry: [6957](#)

UniProt Entry: [P0DSE2](#)

Application Information

Clonality: Armenian Hamster monoclonal antibody

Clone ID: 24GB15645

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: FITC-Cy5

Immunogen

The immunogen information is proprietary.

Isotype

Armenian Hamster IgG

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 μ L/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-TCR β Armenian Hamster Monoclonal Antibody, FITC-Cy5 Conjugated

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM