

Anti-FOXP3 Rat Monoclonal Antibody, AF647 Conjugated



Catalog #: F0786-AF647

Aliases

FOXP3; Forkhead Box P3; DIETER; XPID; AIID; PIDX; IPEX; JM2; Immune Dysregulation, Polyendocrinopathy, Enteropathy, X-Linked; Forkhead Box Protein P3; SCURFIN; Immunodeficiency, Polyendocrinopathy, Enteropathy, X-Linked; FOXP3delta7; Scurfin

Background

Gene Name: FOXP3

NCBI Gene Entry: [50943](#)

UniProt Entry: [Q9BZS1](#)

Application Information

Clonality: Rat monoclonal antibody

Clone ID: 24GB15610

Species Reactivity: Mouse

Applications Tested: Intracellular Staining for Flow Cytometry ☐ ICFC ☐

Conjugation: AF647

Immunogen

The immunogen information is proprietary.

Isotype

Rat IgG2a

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

