

Anti-CD49b Rat Monoclonal Antibody, APC-Cy7 Conjugated



Catalog #: F0794-AC7

Aliases

ITGA2; Integrin Subunit Alpha 2; GPIa; VLAA2; HPA-5; CD49B; Integrin, Alpha 2 (CD49B, Alpha 2 Subunit Of VLA-2 Receptor); Very Late Activation Receptor Subunit Alpha-2; Alpha 2 Subunit Of VLA-2 Receptor; CD49 Antigen-Like Family Member B; Platelet Membrane Glycoprotein Ia; Human Platelet Alloantigen 5; Collagen Receptor; Alpha 2 Integrin; Integrin Alpha-2; Glycoprotein Ia; Very Late Activation Protein 2 Receptor, Alpha-2 Subunit; Very Late Activation Receptor Alpha-2 Subunit; Human Platelet Alloantigen System 5; Platelet Glycoprotein GPIa; Platelet Antigen Br; VLA-2 Subunit Alpha; CD49b Antigen; VLA-2; BR

Background

Gene Name: ITGA2

NCBI Gene Entry: [3673](#)

UniProt Entry: [P17301](#)

Application Information

Clonality: Rat monoclonal antibody

Clone ID: 24GB15585

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: APC-Cy7

Immunogen

The immunogen information is proprietary.

Isotype

Rat IgM

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041


WWW.GENUINBIOTECH.COM

Anti-CD49b Rat Monoclonal Antibody, APC-Cy7 Conjugated

PAGE 2

Note: This product is for research use only.

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM