

Anti-CD154 Armenian Hamster Monoclonal Antibody, PE Conjugated



Catalog #: F0800-PE

Aliases

CD40LG; CD40 Ligand; CD40-L; CD40L; TRAP; TNF-Related Activation Protein; HCD40L; TNFSF5; CD154; HIGM1; T-BAM; Gp39; Tumor Necrosis Factor (Ligand) Superfamily Member 5; T-B Cell-Activating Molecule; CD40 Antigen Ligand; T-Cell Antigen Gp39; IMD3; Tumor Necrosis Factor (Ligand) Superfamily, Member 5 (Hyper-IgM Syndrome); Tumor Necrosis Factor Ligand Superfamily Member 5; Hyper-IgM Syndrome; CD154 Antigen; IGM

Background

Gene Name: CD40LG

NCBI Gene Entry: [959](#)

UniProt Entry: [P29965](#)

Application Information

Clonality: Armenian Hamster monoclonal antibody

Clone ID: 24GB15670

Species Reactivity: Mouse

Applications Tested: Flow Cytometry (FCM)

Conjugation: PE

Immunogen

The immunogen information is proprietary.

Isotype

Armenian Hamster IgG3

Storage Buffer

Supplied in PBS (pH 7.2-7.4) containing 0.2% BSA and 0.1% Proclin 300.

Storage

Store at 4 °C and protected from prolonged exposure to light.

Recommended Dilutions

5 µL/Tests

Note: This product is for research use only.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS


SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-CD154 Armenian Hamster Monoclonal Antibody, PE Conjugated

PAGE 2

Validation Data



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM