

Anti-Phospho-c-CBL (S669) Rabbit Monoclonal Antibody



Catalog #: U0103

Aliases

CBL2; RNF55; E3 ubiquitin-protein ligase CBL; Casitas B-lineage lymphoma proto-oncogene; Proto-oncogene c-Cbl; RING finger protein 55; Signal transduction protein CBL

Background

Gene Name: CBL

NCBI Gene Entry: [867](#)

UniProt Entry: [P22681](#)

Application Information

Molecular Weight: Predicted, 99 kDa; observed, 120 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7520

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-c-CBL (S669)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of c-CBL (Phospho-S669) expression in HeLa (A) whole cell lysates.