

Anti-Phospho-Nucleolin (T76) Rabbit Monoclonal Antibody



Catalog #: U0105

Aliases

Nucleolin; Protein C23

Background

Gene Name: NCL

NCBI Gene Entry: [4691](#)

UniProt Entry: [P19338](#)

Application Information

Molecular Weight: Predicted, 76 kDa; observed, 100 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7530

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Nucleolin (T76)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-Nucleolin (T76) Rabbit Monoclonal Antibody

PAGE 2



Western blotting analysis of Nucleolin (Phospho-T76) expression in K562 (A) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM