

Anti-Phospho-c-Myc (T58+S62) Rabbit Monoclonal Antibody



Catalog #: U0106

Aliases

BHLHE39; Myc proto-oncogene protein; Class E basic helix-loop-helix protein 39; bHLHe39; Proto-oncogene c-Myc; Transcription factor p64

Background

Gene Name: MYC

NCBI Gene Entry: [4609](#)

UniProt Entry: [P01106](#)

Application Information

Molecular Weight: Predicted, 48 kDa; observed, 57 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7535

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-c-Myc (T58+S62)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-c-Myc (T58+S62) Rabbit Monoclonal Antibody

PAGE 2



Western blotting analysis of c-Myc (Phospho-T58+S62) expression in THP1 (A), HeLa (B), mouse lung (C) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM