

Anti-Phospho-PRKAR2A (S99) Rabbit Monoclonal Antibody



Catalog #: U0107

Aliases

PKR2; PRKAR2; cAMP-dependent protein kinase type II-alpha regulatory subunit

Background

Gene Name: PRKAR2A

NCBI Gene Entry: [5576](#)

UniProt Entry: [P13861](#)

Application Information

Molecular Weight: Predicted, 45 kDa; observed, 51 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7540

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-PRKAR2A (S99)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-PRKAR2A (S99) Rabbit Monoclonal Antibody

PAGE 2



Western blotting analysis of PRKAR2A (Phospho-S99) expression in mouse liver (A), rat liver (B) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM