

# Anti-Phospho-PPP2CA (Y307) Rabbit Monoclonal Antibody



**Catalog #: U0108**

## Aliases

Serine/threonine-protein phosphatase 2A catalytic subunit alpha isoform; PP2A-alpha; Replication protein C; RP-C

## Background

Gene Name: PPP2CA

NCBI Gene Entry: [5515](#)

UniProt Entry: [P67775](#)

## Application Information

Molecular Weight: Predicted, 35 kDa; observed, 35 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7545

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

## Immunogen

A synthesized peptide derived from human Phospho-PPP2CA (Y307)

## Isotype

Rabbit IgG

## Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

**Note:** This product is for research use only.

## Validation Data

---


### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Western blotting analysis of PPP2CA (Phospho-Y307) expression in HEK293T (A), Hela (B), A431 (C), mouse muscle (D) whole cell lysates.