

Anti-Phospho-Histone H3 (S31) Rabbit Monoclonal Antibody



Catalog #: U0110

Aliases

H3.3A; H3F3; H3.3B; Histone H3.3

Background

Gene Name: H3F3A/B

NCBI Gene Entry: [3020/3021](#)

UniProt Entry: [P84243](#)

Application Information

Molecular Weight: Predicted, 15 kDa; observed, 15 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7555

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Histone H3 (S31)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-Histone H3 (S31) Rabbit Monoclonal Antibody

PAGE 2



Western blotting analysis of Histone H3 (Phospho-S31) expression in HeLa (A), HuT78 (B), MCF7 (C) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM