

Anti-Phospho-FOXO3 (S253) Rabbit Monoclonal Antibody



Catalog #: U0115

Aliases

FKHRL1; FOXO3A; Forkhead box protein O3; AF6q21 protein; Forkhead in rhabdomyosarcoma-like 1

Background

Gene Name: FOXO3

NCBI Gene Entry: [2309](#)

UniProt Entry: [O43524](#)

Application Information

Molecular Weight: Predicted, 71 kDa; observed, 97 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7580

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-FOXO3 (S253)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of FOXO3 (Phospho-S253) expression in HEK293T (A), SHSY5Y (B), A431 (C), mouse liver (D), mouse kidney (E), rat liver (F) whole cell lysates.