

Anti-Phospho-STAT5 (Y694) Rabbit Monoclonal Antibody



Catalog #: U0116

Aliases

STAT5A; STAT5; Signal transducer and activator of transcription 5A; STAT5B; Signal transducer and activator of transcription 5B

Background

Gene Name: STAT5A/B

NCBI Gene Entry: [6776/6777](#)

UniProt Entry: [P42229/P51692](#)

Application Information

Molecular Weight: Predicted, 90,89 kDa; observed, 90 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7585

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-STAT5 (Y694)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-STAT5 (Y694) Rabbit Monoclonal Antibody

PAGE 2



Western blotting analysis of STAT5 (Phospho-Y694) expression in K562 (A), mouse liver (B) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM