

Anti-Phospho-IRE1 (S724) Rabbit Monoclonal Antibody



Catalog #: U0118

Aliases

IRE1; Serine/threonine-protein kinase/endoribonuclease IRE1; Endoplasmic reticulum-to-nucleus signaling 1; Inositol-requiring protein 1; hIRE1p; Ire1-alpha; IRE1a

Background

Gene Name: ERN1

NCBI Gene Entry: [2081](#)

UniProt Entry: [O75460](#)

Application Information

Molecular Weight: Predicted, 109 kDa; observed, 110 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7595

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-IRE1 (S724)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of IRE1 (Phospho-S724) expression in K562 (A) whole cell lysates.