

Anti-Phospho-Synapsin 1 (S9) Rabbit Monoclonal Antibody



Catalog #: U0119

Aliases

Synapsin-1; Brain protein 4.1; Synapsin I

Background

Gene Name: SYN1

NCBI Gene Entry: [6853](#)

UniProt Entry: [P17600](#)

Application Information

Molecular Weight: Predicted, 74 kDa; observed, 77 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7600

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Synapsin 1 (S9)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-Synapsin 1 (S9) Rabbit Monoclonal Antibody

PAGE 2



Western blotting analysis of Synapsin 1 (Phospho-S9) expression in mouse brain (A) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM