

Anti-Phospho-RPA2 (T21) Rabbit Monoclonal Antibody



Catalog #: U0121

Aliases

REPA2; RPA32; RPA34; Replication protein A 32 kDa subunit; RP-A p32; Replication factor A protein 2; RF-A protein 2; Replication protein A 34 kDa subunit; RP-A p34

Background

Gene Name: RPA2

NCBI Gene Entry: [6118](#)

UniProt Entry: [P15927](#)

Application Information

Molecular Weight: Predicted, 29 kDa; observed, 32 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7610

Species Reactivity: Human, mouse

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-RPA2 (T21)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of RPA2 (Phospho-T21) expression in k562 (A), HeLa (B), PC3 (C), mouse kidney (D), mouse muscle (E) whole cell lysates.