

Anti-Phospho-SMC1A (S957) Rabbit Monoclonal Antibody



Catalog #: U0122

Aliases

DXS423E; KIAA0178; SB1.8; SMC1; SMC1L1; Structural maintenance of chromosomes protein 1A; SMC protein 1A; SMC-1-alpha; SMC-1A; Sb1.8

Background

Gene Name: SMC1A

NCBI Gene Entry: [8243](#)

UniProt Entry: [Q14683](#)

Application Information

Molecular Weight: Predicted, 143 kDa; observed, 160 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7615

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-SMC1A (S957)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of SMC1A (Phospho-S957) expression in MCF7 (A), A431 (B) whole cell lysates.