

Anti-Phospho-HSL (S853) Rabbit Monoclonal Antibody



Catalog #: U0123

Aliases

Hormone-sensitive lipase; HSL

Background

Gene Name: LIPE

NCBI Gene Entry: [3991](#)

UniProt Entry: [Q05469](#)

Application Information

Molecular Weight: Predicted, 116 kDa; observed, 83 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7620

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-HSL (S853)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of HSL (Phospho-S853) expression in mouse muscle (A), rat muscle (B) whole cell lysates.