

Anti-Phospho-MKK1 (T292) Rabbit Monoclonal Antibody



Catalog #: U0127

Aliases

MEK1; PRKMK1; Dual specificity mitogen-activated protein kinase kinase 1; MAP kinase kinase 1; MAPKK 1; MKK1; ERK activator kinase 1; MAPK/ERK kinase 1; MEK 1

Background

Gene Name: MAP2K1

NCBI Gene Entry: [5604](#)

UniProt Entry: [Q02750](#)

Application Information

Molecular Weight: Predicted, 43 kDa; observed, 45 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7640

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-MKK1 (T292)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of MKK1 (Phospho-T292) expression in A431 (A), A549 (B), mouse muscle (C), rat muscle (D) whole cell lysates.