

Anti-Phospho-PAK4/5/6 (S474/560/602) Rabbit Monoclonal Antibody



Catalog #: U0137

Aliases

KIAA1142; Serine/threonine-protein kinase PAK 4; p21-activated kinase 4; PAK-4; PAK5; Serine/threonine-protein kinase PAK 6; PAK-5; p21-activated kinase 6; PAK-6

Background

Gene Name: PAK4/5/6

NCBI Gene Entry: [10298/56924](#)

UniProt Entry: [O96013/Q9P286/Q9NQU5](#)

Application Information

Molecular Weight: Predicted, 64 kDa; observed, 72 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7690

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-PAK4/5/6 (S474/560/602)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-PAK4/5/6 (S474/560/602) Rabbit Monoclonal Antibody

PAGE 2



Western blotting analysis of PAK4/5/6 (Phospho-S474/560/602) expression in K562 (A), MCF7 (B), PC3 (C) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM