

Anti-Phospho-TAK1 (S439) Rabbit Monoclonal Antibody



Catalog #: U0141

Aliases

TAK1; Mitogen-activated protein kinase kinase kinase 7; Transforming growth factor-beta-activated kinase 1; TGF-beta-activated kinase 1

Background

Gene Name: MAP3K7

NCBI Gene Entry: [6885](#)

UniProt Entry: [O43318](#)

Application Information

Molecular Weight: Predicted, 67 kDa; observed, 72 kDa

Clonality: Rabbit monoclonal antibody

Clone ID: 24GB7710

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-TAK1 (S439)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of TAK1 (Phospho-S439) expression in PC3 (A), HeLa (B), K562 (C) whole cell lysates.