

# Anti-Phospho-CDC25A (S124) Rabbit Polyclonal Antibody



**Catalog #: U0166**

## Aliases

M-phase inducer phosphatase 1; Dual specificity phosphatase Cdc25A

## Background

Gene Name: CDC25A

NCBI Gene Entry: [993](#)

UniProt Entry: [P30304](#)

## Application Information

Molecular Weight: Predicted, 59 kDa; observed, 70 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine

Applications Tested: Western blotting (WB)

## Immunogen

A synthesized peptide derived from human Phospho-CDC25A (S124)

## Isotype

Rabbit IgG

## Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

**Note:** This product is for research use only.

## Validation Data

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041


### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)

## Anti-Phospho-CDC25A (S124) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of CDC25A (Phospho-S124) expression in K562 (A), HEK293T (B) whole cell lysates.

---

### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)