

Anti-Phospho-C/EBP alpha (T226) Rabbit Polyclonal Antibody



Catalog #: U0171

Aliases

CCAAT/enhancer-binding protein alpha; C/EBP alpha

Background

Gene Name: CEBPA

NCBI Gene Entry: [1050](#)

UniProt Entry: [P49715](#)

Application Information

Molecular Weight: Predicted, 37 kDa; observed, 45 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, chicken, zebrafish

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-C/EBP alpha (T226)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-C/EBP alpha (T226) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of C/EBP alpha (Phospho-T226) expression in MCF7 EGF-treated (A), mouse spleen (B), mouse liver (C), rat spleen (D), rat liver (E) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM