

Anti-Phospho-E2F1 (T433) Rabbit Polyclonal Antibody



Catalog #: U0185

Aliases

RBBP3; Transcription factor E2F1; E2F-1; PBR3; Retinoblastoma-associated protein 1; RBAP-1; Retinoblastoma-binding protein 3; RBBP-3; pRB-binding protein E2F-1

Background

Gene Name: E2F1

NCBI Gene Entry: [1869](#)

UniProt Entry: [Q01094](#)

Application Information

Molecular Weight: Predicted, 46 kDa; observed, 47 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-E2F1 (T433)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of E2F1 (Phospho-T433) expression in Jurkat colchicine-treated (A), mouse spleen (B), PC12 colchicine-treated (C) whole cell lysates.



Lane A,B,C,D : AML12 cells treated with PMA (0.1 nmol)

Western blotting analysis of E2F1 (Phospho-T433) expression in AML12-PMA treated (A) (B) (C) (D) whole cell lysates.