

Anti-Phospho-GAB1 (Y659) Rabbit Polyclonal Antibody



Catalog #: U0201

Aliases

GRB2-associated-binding protein 1; GRB2-associated binder 1; Growth factor receptor bound protein 2-associated protein 1

Background

Gene Name: GAB1

NCBI Gene Entry: [2549](#)

UniProt Entry: [Q13480](#)

Application Information

Molecular Weight: Predicted, 76 kDa; observed, 110 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-GAB1 (Y659)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of GAB1 (Phospho-Y659) expression in HeLa insulin-treated (A), NIH3T3 (B), C6 (C), Raw264.7 insulin-treated (D), H9C2 insulin-treated (E) whole cell lysates.