

Anti-Phospho-GRLF1 (Y1087) Rabbit Polyclonal Antibody



Catalog #: U0205

Aliases

GRF1; GRLF1; KIAA1722; Rho GTPase-activating protein 35; Glucocorticoid receptor DNA-binding factor 1; Glucocorticoid receptor repression factor 1; GRF-1; Rho GAP p190A; p190-A

Background

Gene Name: ARHGAP35

NCBI Gene Entry: [2909](#)

UniProt Entry: [Q9NRY4](#)

Application Information

Molecular Weight: Predicted, 170 kDa; observed, 190 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, dog, pig

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-GRLF1 (Y1087)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of GRLF1 (Phospho-Y1087) expression in mouse testis (A), MCF7 (B), HUT78 (C), K562 (D) whole cell lysates.