

Anti-Phospho-CD54 (Y512) Rabbit Polyclonal Antibody



Catalog #: U0211

Aliases

Intercellular adhesion molecule 1; ICAM-1; Major group rhinovirus receptor; CD54

Background

Gene Name: ICAM1

NCBI Gene Entry: [3383](#)

UniProt Entry: [P05362](#)

Application Information

Molecular Weight: Predicted, 57 kDa; observed, 92 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-CD54 (Y512)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of CD54 (Phospho-Y512) expression in Saos2 (A), HepG2 TNFa-treated (B), Raji (C) whole cell lysates.