

Anti-Phospho-MARCKS (S163) Rabbit Polyclonal Antibody



Catalog #: U0233

Aliases

MACS; PRKCSL; Myristoylated alanine-rich C-kinase substrate; MARCKS; Protein kinase C substrate, 80 kDa protein, light chain; 80K-L protein; PKCSL

Background

Gene Name: MARCKS

NCBI Gene Entry: [4082](#)

UniProt Entry: [P29966](#)

Application Information

Molecular Weight: Predicted, 31 kDa; observed, 75 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, chicken, zebrafish

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-MARCKS (S163)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of MARCKS (Phospho-S163) expression in HeLa (A), Lovo (B), HepG2 (C) whole cell lysates.