

Anti-Phospho-FOXO4 (T451) Rabbit Polyclonal Antibody



Catalog #: U0242

Aliases

AFX; AFX1; MLLT7; Forkhead box protein O4; Fork head domain transcription factor AFX1

Background

Gene Name: FOXO4

NCBI Gene Entry: [4303](#)

UniProt Entry: [P98177](#)

Application Information

Molecular Weight: Predicted, 53 kDa; observed, 65 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-FOXO4 (T451)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of FOXO4 (Phospho-T451) expression in THP1 EGF-treated (A), K562 (B), Jurkat (C), NIH3T3 (D), mouse lung (E) whole cell lysates.