

Anti-Phospho-Nucleophosmin (T234) Rabbit Polyclonal Antibody



Catalog #: U0254

Aliases

NPM; Nucleophosmin; NPM; Nucleolar phosphoprotein B23; Nucleolar protein NO38; Numatrin

Background

Gene Name: NPM1

NCBI Gene Entry: [4869](#)

UniProt Entry: [P06748](#)

Application Information

Molecular Weight: Predicted, 32 kDa; observed, 38 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, rabbit

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Nucleophosmin (T234)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Nucleophosmin (Phospho-T234) expression in MCF7 colchicine-treated (A), LNCaP (B), HeLa (C), Raw264.7 colchicine-treated (D), PC12 colchicine-treated (E) whole cell lysates.