

Anti-Phospho-PPAR gamma (S112) Rabbit Polyclonal Antibody



Catalog #: U0269

Aliases

NR1C3; Peroxisome proliferator-activated receptor gamma; PPAR-gamma; Nuclear receptor subfamily 1 group C member 3

Background

Gene Name: PPARG

NCBI Gene Entry: [5468](#)

UniProt Entry: [P37231](#)

Application Information

Molecular Weight: Predicted, 57 kDa; observed, 54, 57 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, dog, rabbit

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-PPAR gamma (S112)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of PPAR gamma (Phospho-S112) expression in mouse heart (A), rat heart (B), H9C2 (C), AML12 (D), A549 (E), A2780 (F), LO2 (G) whole cell lysates.