

Anti-Phospho-PKC delta (Y52) Rabbit Polyclonal Antibody



Catalog #: U0273

Aliases

Protein kinase C delta type; Tyrosine-protein kinase PRKCD; nPKC-delta

Background

Gene Name: PRKCD

NCBI Gene Entry: [5580](#)

UniProt Entry: [Q05655](#)

Application Information

Molecular Weight: Predicted, 77 kDa; observed, 78 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, chicken, dog, monkey, pig, rabbit

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-PKC delta (Y52)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of PKC delta (Phospho-Y52) expression in MCF7 serum starvation-treated (A), HeLa PMA-treated (B), Jurkat (C), Raw264.7 serum starvation-treated (D), H9C2 serum starvation-treated (E) whole cell lysates.