

Anti-Phospho-PKC mu (S910) Rabbit Polyclonal Antibody



Catalog #: U0274

Aliases

PKD; PKD1; PRKCM; Serine/threonine-protein kinase D1; Protein kinase C mu type; Protein kinase D; nPKC-D1; nPKC-mu

Background

Gene Name: PRKD1

NCBI Gene Entry: [5587](#)

UniProt Entry: [Q15139](#)

Application Information

Molecular Weight: Predicted, 101 kDa; observed, 115 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-PKC mu (S910)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-PKC mu (S910) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of PKC mu (Phospho-S910) expression in HeLa (A), Jurkat (B), HEK293T PMA-treated (C), NIH3T3 PDGF-treated (D), H9C2 PMA-treated (E) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM