

Anti-Phospho-Focal Adhesion Kinase (Y576) Rabbit Polyclonal Antibody



Catalog #: U0280

Aliases

FAK; FAK1; Focal adhesion kinase 1; FADK 1; Focal adhesion kinase-related nonkinase; FRNK; Protein phosphatase 1 regulatory subunit 71; PPP1R71; Protein-tyrosine kinase 2; p125FAK; pp125FAK

Background

Gene Name: PTK2

NCBI Gene Entry: [5747](#)

UniProt Entry: [Q05397](#)

Application Information

Molecular Weight: Predicted, 119, 99, 63, 48, 120, 39, 114 kDa; observed, 125 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, chicken

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Focal Adhesion Kinase (Y576)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-Focal Adhesion Kinase (Y576) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of Focal Adhesion Kinase (Phospho-Y576) expression in HEK293T Forskolin-treated (A), NIH3T3 Forskolin-treated (B), H9C2 Forskolin-treated (C) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM