

Anti-Phospho-SHB (Y246) Rabbit Polyclonal Antibody



Catalog #: U0294

Aliases

SH2 domain-containing adapter protein B

Background

Gene Name: SHB

NCBI Gene Entry: [6461](#)

UniProt Entry: [Q15464](#)

Application Information

Molecular Weight: Predicted, 55 kDa; observed, 55 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-SHB (Y246)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of SHB (Phospho-Y246) expression in Hela (A), A375 (B), rat kidney (C), rat spleen (D) whole cell lysates.