

Anti-Phospho-TERF1 (S219) Rabbit Polyclonal Antibody



Catalog #: U0307

Aliases

PIN2; TRBF1; TRF; TRF1; Telomeric repeat-binding factor 1; NIMA-interacting protein 2; TTAGGG repeat-binding factor 1; Telomeric protein Pin2/TRF1

Background

Gene Name: TERF1

NCBI Gene Entry: [7013](#)

UniProt Entry: [P54274](#)

Application Information

Molecular Weight: Predicted, 50 kDa; observed, 60 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-TERF1 (S219)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of TERF1 (Phospho-S219) expression in DLD UV-treated (A), Raw264.7 UV-treated (B), PC12 UV-treated (C) whole cell lysates.