

Anti-Phospho-Tensin 1 (Y1326) Rabbit Polyclonal Antibody



Catalog #: U0308

Aliases

TNS; Tensin-1

Background

Gene Name: TNS1

NCBI Gene Entry: [7145](#)

UniProt Entry: [Q9HBL0](#)

Application Information

Molecular Weight: Predicted, 185 kDa; observed, 200 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, pig

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Tensin 1 (Y1326)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Tensin 1 (Phospho-Y1326) expression in MEF (A), PC12 (B), LOVO (C) whole cell lysates.