

Anti-Phospho-MAX (S11) Rabbit Polyclonal Antibody



Catalog #: U0340

Aliases

BHLHD4; Protein max; Class D basic helix-loop-helix protein 4; bHLHd4; Myc-associated factor X

Background

Gene Name: MAX

NCBI Gene Entry: [4149](#)

UniProt Entry: [P61244](#)

Application Information

Molecular Weight: Predicted, 18 kDa; observed, 18 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, chicken

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-MAX (S11)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of MAX (Phospho-S11) expression in Hela (A), mouse heart (B), H9C2 (C) whole cell lysates.