

Anti-Phospho-ROS (Y2114) Rabbit Polyclonal Antibody



Catalog #: U0346

Aliases

MCF3; ROS; Proto-oncogene tyrosine-protein kinase ROS; Proto-oncogene c-Ros; Proto-oncogene c-Ros-1; Receptor tyrosine kinase c-ros oncogene 1; c-Ros receptor tyrosine kinase

Background

Gene Name: ROS1

NCBI Gene Entry: [6098](#)

UniProt Entry: [P08922](#)

Application Information

Molecular Weight: Predicted, 263 kDa; observed, 263 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, bovine, chicken, pig, rabbit

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-ROS (Y2114)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ROS (Phospho-Y2114) expression in SKOV3A (A), Raw264.7 (B), PC12 (C) whole cell lysates.