

# Anti-Phospho-TLK1 (S743) Rabbit Polyclonal Antibody



**Catalog #: U0347**

## Aliases

TLK1; KIAA0137; Serine/threonine-protein kinase tousled-like 1; PKU-beta; Tousled-like kinase 1

## Background

Gene Name: TLK1

NCBI Gene Entry: [9874](#)

UniProt Entry: [Q9UKI8](#)

## Application Information

Molecular Weight: Predicted, 86 kDa; observed, 86 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, rat, monkey

Applications Tested: Western blotting (WB)

## Immunogen

A synthesized peptide derived from human Phospho-TLK1 (S743)

## Isotype

Rabbit IgG

## Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

## Storage

Store at -20 °C for one year.

## Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

**Note:** This product is for research use only.

## Validation Data

---


### SUPPORT

SUPPORT@GENUINBIOTECH.COM  
TEL: +1-540-855-7041

### ORDERS

SALES@GENUINBIOTECH.COM  
FAX: +1-540-855-7041

[WWW.GENUINBIOTECH.COM](http://WWW.GENUINBIOTECH.COM)



Western blotting analysis of TLK1 (Phospho-S743) expression in MCF7 (A), rat brain (B), H9C2 (C) whole cell lysates.