

Anti-Phospho-HSL (S855) Rabbit Polyclonal Antibody



Catalog #: U0385

Aliases

Hormone-sensitive lipase; HSL

Background

Gene Name: LIPE

NCBI Gene Entry: [3991](#)

UniProt Entry: [Q05469](#)

Application Information

Molecular Weight: Predicted, 116 kDa; observed, 81, 83 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, pig

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-HSL (S855)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of HSL (Phospho-S855) expression in Human muscle (A) whole cell lysates.