

Anti-Phospho-Nibrin (S343) Rabbit Polyclonal Antibody



Catalog #: U0395

Aliases

NBS; NBS1; P95; Nibrin; Cell cycle regulatory protein p95; Nijmegen breakage syndrome protein 1

Background

Gene Name: NBN

NCBI Gene Entry: [4683](#)

UniProt Entry: [O60934](#)

Application Information

Molecular Weight: Predicted, 84 kDa; observed, 95 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-Nibrin (S343)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of Nibrin (Phospho-S343) expression in Jurkat (A), LO2 UV-treated (B), Raw264.7 UV-treated (C), PC12 (D) whole cell lysates.