

Anti-Phospho-EPHA7 (Y608) Rabbit Polyclonal Antibody



Catalog #: U0466

Aliases

EHK3; HEK11; Ephrin type-A receptor 7; EPH homology kinase 3; EHK-3; EPH-like kinase 11; EK11; hEK11

Background

Gene Name: EPHA7

NCBI Gene Entry: [2045](#)

UniProt Entry: [Q15375](#)

Application Information

Molecular Weight: Predicted, 112 kDa; observed, 130 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, chicken

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-EPHA7 (Y608)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of EPHA7 (Phospho-Y608) expression in A549 (A), NS-1 (B), PC12 (C) whole cell lysates.