

Anti-Phospho-ATM (S1981) Rabbit Polyclonal Antibody



Catalog #: U0539

Aliases

Serine-protein kinase ATM; Ataxia telangiectasia mutated; A-T mutated

Background

Gene Name: ATM

NCBI Gene Entry: [472](#)

UniProt Entry: [Q13315](#)

Application Information

Molecular Weight: Predicted, 350 kDa; observed, 350 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-ATM (S1981)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of ATM (Phospho-S1981) expression in K562 (A), HEK293T UV-treated (B) whole cell lysates.