

Anti-Phospho-BCL6 (S333) Rabbit Polyclonal Antibody



Catalog #: U0544

Aliases

BCL5; LAZ3; ZBTB27; ZNF51; B-cell lymphoma 6 protein; BCL-6; B-cell lymphoma 5 protein; BCL-5; Protein LAZ-3; Zinc finger and BTB domain-containing protein 27; Zinc finger protein 51

Background

Gene Name: BCL6

NCBI Gene Entry: [604](#)

UniProt Entry: [P41182](#)

Application Information

Molecular Weight: Predicted, 78 kDa; observed, 79 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, dog

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-BCL6 (S333)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data


SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM



Western blotting analysis of BCL6 (Phospho-S333) expression in Jurkat (A), mouse brain (B) whole cell lysates.