

Anti-Phospho-CDC25B (S323) Rabbit Polyclonal Antibody



Catalog #: U0549

Aliases

CDC25HU2; M-phase inducer phosphatase 2; Dual specificity phosphatase Cdc25B

Background

Gene Name: CDC25B

NCBI Gene Entry: [994](#)

UniProt Entry: [P30305](#)

Application Information

Molecular Weight: Predicted, 64 kDa; observed, 65 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, monkey

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-CDC25B (S323)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-CDC25B (S323) Rabbit Polyclonal Antibody

PAGE 2



Western blotting analysis of CDC25B (Phospho-S323) expression in A549 (A), NIH3T3 (B) whole cell lysates.

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM