

Anti-Phospho-PTP alpha (Y798) Rabbit Polyclonal Antibody



Catalog #: U0581

Aliases

PTPA; PTPRL2; Receptor-type tyrosine-protein phosphatase alpha; Protein-tyrosine phosphatase alpha; R-PTP-alpha

Background

Gene Name: PTPRA

NCBI Gene Entry: [5786](#)

UniProt Entry: [P18433](#)

Application Information

Molecular Weight: Predicted, 90 kDa; observed, 120 kDa

Clonality: Rabbit polyclonal antibody

Species Reactivity: Human, mouse, rat, chicken, dog, monkey, rabbit, sheep, zebrafish

Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human Phospho-PTP alpha (Y798)

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.3) containing 30% glycerol, and 0.01% sodium azide.

Storage

Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:500-1:1,000

Note: This product is for research use only.

Validation Data

SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041


ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM

Anti-Phospho-PTP alpha (Y798) Rabbit Polyclonal Antibody

PAGE 2



SUPPORT

SUPPORT@GENUINBIOTECH.COM
TEL: +1-540-855-7041

ORDERS

SALES@GENUINBIOTECH.COM
FAX: +1-540-855-7041

WWW.GENUINBIOTECH.COM